CalFast

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Bibliografia

1. Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
2. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol. 1992;27:793-798.
3. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
4. Tibble J. et al.: A simple method for assessing intestinal

- inflammation in Crohn's disease. Gut 2000;47:506-513.
5. Campeotto F. et al.: Fecal calprotectin: cutoff values for identifying intestinal distress in preterm infants. J Pediatr Gastroenterol Nutr. 2009 Apr;48(4):507-10.
6. Sohnle PG. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.
7. Steinbakk M. Et al.: Antimicrobial actions of calcium binding leukocyte L1 protein, calprotectin. Lancet 1990;336:763-765.
8. Ton H. et al.: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000;292:41-54.
9. Røseth A.G. et al.: Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.
10. Johnne B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

CalFast

20 test, codice 9191

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2013.11.11

Eurospital[®]

Eurospital SpA

Via Flavia 122, 34147 Trieste, Italia

Tel. +39 040 8997.1 Fax +39 040 280944

www.eurospital.com - info@eurospital.com

CalFast

Rapid test for the quantitative detection of calprotectin in stools.

For in vitro diagnostic use.

Intended use

CalFast is an immunochromatographic assay for the quantitative detection of calprotectin in faeces. The calprotectin extraction occurs in a dedicated device, which allows to greatly reduce the time needed for the pre-analytical procedure.

Test Principle

CalFast is immunochromatographic assay which combines a mixture of anti-calprotectin polyclonal and monoclonal antibodies. Calprotectin present in the diluted sample is bound by the gold-conjugated antibody. The so formed complex moves along the reaction strip and binds to the anti-calprotectin antibodies present in the reaction area (T) while unbound antibodies will bind in the control area (C). The concentration of calprotectin in the sample is calculated by means of the CalFast reader.

Provided materials and preparation

Cassettes	20 pieces
Diluent	1 vial
Stool collection and extraction devices	20 pieces

Composition of provided materials/reagents

Cassettes

Twenty (20) Cassettes single-packed in a sealed aluminum foil containing a desiccant. Cassettes are disposable.

Diluent

One (1) Vial containing 50 ml. Ready to use.

Stool collection and extraction devices

Twenty (20) Devices containing a ready-to-use extraction solution. Devices are disposable.

Materials required but not provided

1. Plastic tubes for faecal extract dilution.
2. Vortex.
3. Precision pipettes with disposable tips (10-200 μ l, 1ml).
4. Centrifuge for microtubes.

5. Eppendorf tubes (or similar) for extract dilution.

6. Timer.

7. CalFast Reader.

Performance criteria and limits of the test

Intra-assay variation

Four (4) stool samples were tested twelve (12) times in the same run. The obtained results are as follows:

Sample	1	2	3	4
Mean mg/Kg	208,8	142,4	101,2	72,2
STD	25,5	15,6	8,5	8,2
CV%	12,2	10,9	8,4	11,3

Inter-assay variation

Five (5) stool samples were tested in duplicate in three (3) different series. The obtained results are as follows:

Sample	1	2	3	4	5
Mean mg/Kg	88,4	98,3	114,4	172,3	200,9
STD	7,6	9,3	7,9	16,7	23,6
CV%	8,7	9,4	6,9	9,7	11,8

Inter-lot variation

Five (5) stool samples were tested in duplicate using three (3) different batches of product. The obtained results are as follows:

Sample	1	2	3	4	5
Mean mg/Kg	84,4	101,0	121,9	104,6	221,7
STD	7,5	8,4	11,5	12,8	22,7
CV%	8,8	8,3	9,4	12,2	10,2

Analytical Sensitivity

A positive samples at a known calprotectin concentration was diluted from 1150 mg/kg to 9 mg/kg and tested accordingly. Each dilution has been tested ten (10) times. Within the assay range (15 -300 mg/kg), all CV% values were lower than the expected 15% value.

Diagnostic sensitivity

Stool samples from twenty (20) subjects affected by clinically diagnosed IBD (Inflammatory Bowel Disease) were tested. The obtained value was 90%.

Diagnostic specificity

Stool samples from hundred and two (102) subjects affected by either IBS (Irritable Bowel Syndrome) or with disease different than IBD were tested. The obtained value was 90,2%.

Concordance between methods

Hundred-sixty-three (163) stool samples from subjects with or without clinical diagnosis of IBD were tested with CalFast (immuno-chromatographic assay). The results, in terms of positive/negative response, were compared with those obtained with Calprest test (ELISA). Hundred-fifty (150) samples showed concordant results (Pos/Pos or Neg/Neg) between the two methods. Thirteen (13) samples showed discordant results (Pos/Neg or Neg/Pos). The overall concordance was 92%.

Storage and expiry date of reagents

All reagents must be stored at 2-8°C and used before expiration date printed on the label.

Stability after first opening

The Sample Diluent is stable for sixty (60) days after first opening, if stored at 2-8°C in its original container.

Transportation stability

An accelerate stability test proved that all components of the kit are stable after storage for 96 hours at 37°C.

Advices and precautions

All patient samples must be treated as potential biohazard, and shall be handled and disposed of according to local laboratory legislation.

Procedure for stool collection

The stool to be tested by CalFast must be collected by means of the Stool Collection Devices included in the kit package.

Sampling procedure for faeces of normal consistency

1. Collect a stool sample in a sterile container.
2. Pull out the shaped stick by turning the screw-cap counter-clockwise.
3. Dip the shaped stick in the previously collected stool sample.
4. Repeatedly rotate the stick in the stool sample until all grooves are filled with faecal material.
5. Before placing the stick into the device, remove the stool excess by rotating the stick on the internal wall of the stool container.
6. Fully insert the stick, with faecal material into the tube containing the extraction solution. Turn the cap clockwise until completely closed.

Sampling procedure for liquid faeces or anyhow having a consistence which does not allow the use of the shake stick

1. It is necessary to have a laboratory pipette.

2. By using the pipette, aspirate 60 µl of liquid stool.
3. Pull out the shaped stick by turning the screw-cap counter-clockwise.
4. Pipette the liquid stool sample into the tube containing the extraction solution.
5. Place the shaped stick back into tube containing the extraction solution. Turn the cap clockwise until completely closed.

Sampling procedure for faeces of hard consistency (the shaped stick can hardly be dipped into stool)

1. Collect a stool sample in a sterile container.
2. Pipette 50-100 µl of saline solution inside the stool container.
3. Let the container with stool and saline solution at room temperature for 60 minutes at least.
4. Proceed then from point 2 of the "Sampling procedure for faeces of normal consistency".

Sampling procedure in case of absence of the device

1. Should the laboratory temporarily not have the collection device, collect a stool sample in a sterile container and keep it at 2-8°C for no more than 7 days. Should the storage period exceed 7 days at 2-8°C, the sample must be frozen at -20°C.
2. Once the device is available again, proceed to the sampling procedure described above. If the sample was frozen, allow it to reach room temperature before to collect the stool sample by means of the device.
3. Proceed from point 2 of the "Sampling procedure for faeces of normal consistency".

Storage of samples

1. Stool can be stored into their original container at 2-8°C up to seven (7) days. For longer storage, it is necessary to store them frozen at -20°C.
2. Once the stool has been collected by means of the stool collection device, the device can be stored at 2-8°C for 48 hours at maximum. After 48 hours, it is mandatory to perform the extraction procedure. On the contrary, the sample must be eliminated.

Storage of the extracts

Once the extraction procedures has occurred, the extracted sample can be immediately used or frozen at -20°C and stored up to three (3) months.

Warnings

1. Wait until all reagents and stool samples reach room temperature (18-28°C) before starting the test.

- In order to obtain reliable quantitative results it is mandatory to centrifuge extracts before proceeding to their dilution.
- Diluted samples must be tested within two (2) hours and cannot be stored for a longer period of time.
- Undiluted extracts removed from the stool collection devices and transferred into clean tubes can be stored for no longer than three (3) months at -20°C.
- It is advised to handle stool as potentially infectious material.

Extraction procedure

- Shake each device by Vortex for 30-60 seconds in order to properly homogenate the content.
- Place the device on a roller shaker for tubes and shake for 20 minutes.
- Repeat steps 1 and 2 for all samples to be processed.
- Remove the devices from the shaker and place them upside-down into in a centrifuge. Centrifuge all devices for 5 minutes at 3000 rpm. Place the device upside down.
- Remove the devices from the centrifuge. Break the low-end and remove it. Transfer the supernatant into a clean tube.

Dilution Procedure


- Dilute the faecal extracts previously prepared 1:50 (e.g. 20µl Sample + 980 µl Diluent Buffer). Mix by vortex (10 sec).
- The Sample is ready to be tested.

Test Procedure

- Open the sealed aluminum foil containing the cassette. Take it and place onto a horizontal surface.
- Take 100 µl of faecal extract and transfer it into the round well of the cassette.
- Incubate at room temperature for 20 minutes.

Reading the results by means of CalFast Reder (same as the picture below)

- Press the "ENTER" key (5) for 1 second to switch on the reader.
- Wait until the reader performs the self-check procedures and calibration control.
- Check that the selected method is <CALFAST>. On the contrary, select <Method> (kind of test) and then <CALFAST>. Press <ENTER>.
- Check the lot number displayed on the reader screen is the same of that report in the CalFast package. On the

contrary, perform the update of the test parameters by means of the "D barcode reader.

- Select <ID Sample> and enter the sample reference number by using the pads on the instrument. Press <ENTER> to confirm.
- Place the processed cassette in the lateral drawer being sure that the round spot is located on the right hand side of the cassette. Close the drawer.
- Select <Measure> and press <ENTER>.
- The CalFast Reader will read the reaction and provide results expressed as mg/Kg.

Warning: For each lot of CalFast, please load the lot specific parameters by reading the 2D barcode located inside the kit package by means of the CalFast barcode reader.

Test parameters are reported in the 2D barcode label printed on the internal part of the cover of the kit: Use of parameters from a different lot number will cause wrong and unreliable results.

Reading of results by means of CalFast Reader New (same as the picture below)

- Check if the CalFast Reader New is connected to a PC or Laptop throughout the dedicated USB cable, included in the box. Turn on the computer, double-click on the icon of CalFast Reader program and enter User and Password.
- Check that all lot info is the same as those reported on the outer label of the CalFast kit. If the lot is different, please update it as described below.
- Place the two-dimensional barcode included in the package into the slot of the instrument.
- Check that the barcode is properly placed inside the Reader correctly.
- Click on this icon  for reading the barcode. All info of the lot will be automatically uploaded into the CalFast Reader New.
- Save the info by clicking on .
- Exit by clicking on this icon .
- This procedure must be repeated every time a new lot is used. All uploaded information for a given lot will be stored until a new batch is used.

Patient's data Input

- Enter at least one of the following data: Family Name, Name, Date of Birth or Patient code.

- Once completed, the green LED lights up, indicating that you can place the cassette inside the slot of the instrument.
- Place the cassette into the slot of the instrument: check if the position is correct by clicking on "Check the strip" (optional).

Check the strip

- Batch code: it indicates the lot currently in use
- Kit name: identification of the product in use.
- Check Strip: to be used before starting a test (optional). When activated, the image of the cassette appears and a check is made on the presence of the Control line. If the Control Line is present and correctly positioned, the screen will display "READABLE". On the contrary, if control line is missing or positioning of the cassette is incorrect, it will show "UNREADABLE".
- Control Strip Picture: it shows the cassette as it is inside the instrument
- Analyzed Image: image of the cassette actually used to perform the test.

Running the test

- If the cassette is readable, start the test by pressing the icon . While reading the cassette, the yellow led turns on and the icon  is replaced by .
- Control Strip box shows two different images: one is the "real" image of the cassette, while the other one is used to perform calculations.
- Analyzed image: it is the image used for calculations. The red and blue lines define the reading range for the calculation.
- CalFast Reader New performs the reading and reports the result in mg/kg or µg/g.
- The concentration value is automatically saved and associated to the specific patient information.
- Warning:** for each lot of CalFast, lot specific parameters must be used. Data can be uploaded by means of the 2D barcode described at step 2. The use of parameters related to a different batch of CalFast can lead to inaccurate results.

For more detailed information, please refer to CalFast Reader New User Manual

Result Validation

- To consider a test result valid, the Control line (C) must be always present. This line is just a functional control of the

test and cannot be used for the interpretation of the Test line (T).

- The amount of calprotectin in the stool sample (T) is calculated by the CalFast Reader.
- If only the line T appears, the test is not valid and must be repeated by using a new cassette.
- If neither the line (C) nor the line (T) show up, the result is not valid and the test must be repeated by using a new cassette.
- NB:** CalFast Reader New can check whether the strip can be read by means of the "Check the Strip".

Results obtained by CalFast Reader

CalFast has been calibrated using both the clinical information of the patients submitted to tests and, at the same time, the results of the same patient samples obtained by Calprest (ELISA).

The CalFast Reader uses a lot specific calibration curve to calculate the Calprotectin concentration. Nonetheless the detectable range of calprotectin is between 15 and 300 mg/kg, it has been decided that the calprotectin values obtained by CalFast reader will be expressed as reported in table below:

Calprotectin values	Calprotectin values on the CalFast reader display
0 - 50 mg/Kg	< 50 mg/Kg
51 - 300 mg/Kg	The real value will be shown
> 300 mg/Kg	> 300 mg/Kg

Reference range and interpretation of results

Good Laboratory Practice suggests that each laboratory should calculate its own normal range. For such calculation, it is requested to test a significant number of samples before assigning a clinical value to this test. An internal analysis provided with the following cut-off values:

Negative: < 70,0 mg/kg
Positive: ≥ 70,0 mg/kg

Samples having a concentration higher than 70 mg/kg must be considered positive to the CalFast test and should be submitted to further clinical investigation to define the reason for the increased Calprotectin concentration.

In healthy adult patients, the average value of Calprotectin is around 25 mg/kg. This value, if obtained by means of CalFast Readers, would be reported as <50 mg/kg. While in adult and pediatric samples, calprotectin values are similar, levels detected in newborns might be significantly higher.

Reference values

Clinical studies (2-5, 10) gave the following values:

Normal healthy reference, median	< 50 mg/kg
Colorectal cancers, median	350 mg/kg
IBD (CD & UC), median	1722 mg/kg

References

1. Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
2. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992;27:793-798.
3. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
4. Tibble J. et al.: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000;47:506-513.
5. Bunn S.K. et al.: Fecal calprotectin: Validation as a noninvasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001;33:14-22.
6. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.
7. Steinbakk M. Et al.: Antimicrobial actions of calcium binding leukocyte L1 protein, calprotectin. Lancet 1990;336:763-765.
8. Ton H. et al.: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000;292:41-54.
9. Røseth A.G. et al.: Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.
10. Johne B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

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References

1. Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
2. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992;27:793-798.
3. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
4. Tibble J. et al.: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000;47:506-513.
5. Bunn S.K. et al.: Fecal calprotectin: Validation as a noninvasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001;33:14-22.
6. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity

of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.

7. Steinbakk M. Et al.: Antimicrobial actions of calcium binding leukocyte L1 protein, calprotectin. Lancet 1990;336:763-765.
8. Ton H. et al.: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000;292:41-54.
9. Røseth A.G. et al.: Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.
10. John B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

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Referencias

1. Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
2. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992;27:793-798.
3. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
4. Tibble J. et al.: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000;47:506-513.
5. Bunn S.K. et al.: Fecal calprotectin: Validation as a noninvasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001;33:14-22.
6. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.
7. Steinbakk M. Et al.: Antimicrobial actions of calcium binding leukocyte L1 protein, calprotectin. Lancet 1990;336:763-765.
8. Ton H. et al.: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000;292:41-54.
9. Røseth A.G. et al.: Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.
10. John B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

Legend / Legenda / Legenda / Légend

CASSETTE	Cassette / Cartucho / Cassetta / Cassette
EXT DEV	Stool Collection and Extraction Device / Dispositivo de recogida y extracción de heces / Dispositivo raccolta ed estrazione feci / Dispositifs de collecte et extraction des selles
DIL	Diluent / Diluyente / Diluente / Diluant

CalFast

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2013.11.11

Eurospital®

Eurospital SpA
Via Flavia 122, 34147 Trieste, Italia
Tel. +39 040 8997.1 Fax +39 040 280944
www.eurospital.com - info@eurospital.com