Calprest®

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Bibliografia

1. Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
2. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992;27:793-798.
3. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
4. Tibble J. et al.: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000;47:506-513.
5. Bunn S.K. et al.: Fecal calprotectin: Validation as a noninvasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001;33:14-22.
6. Dale I. et al.: Purification and partial characterization of a highly immunogenic human leukocyte protein, the L1 antigen. Eur J Biochem 1983;134:1-6.
7. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.
8. Brandtzaeg P. et al.: Distribution of a formalin resistant myelomonocytic antigen (L1) in human tissues. II. Normal and aberrant occurrence in various epithelia. American J of Clin Pathology 1987;87:700-707.
9. Fagerhol M.K.: Nomenclature for proteins: is calprotectin a proper name for the elusive myelomonocytic protein? J Clin Pathol: Mol Pathos 1996;49:M74-M79.

10. Isaksen B. and Fagerhol M.K.: Calprotectin inhibits matrix metalloproteinases by sequestration of zinc. J Clin Pathol: Mol Pathol 2001;54:289-292.
11. Steinbakk M. Et al.: Antimicrobial actions of calcium binding leukocyte L1 protein, calprotectin. Lancet 1990;336:763-765.
12. Yui S. et al.: Induction of apoptotic cell death in mouse lymphoma and human leukemia cell lines by a calcium-binding protein complex, calprotectin, derived from inflammatory peritoneal exudate cells. Journal of Leukocyte Biology; 58:650-658
13. John B. et al.: Functional and clinical aspects of the myelomonocyte protein calprotectin. J Clin Pathol: Molec Pathol 1997;50:113-123.
14. Ton H. et al.: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000;292:41-54.
15. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000;95:2831-2837
16. Røseth A.G. et al.: Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.
17. John B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

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Calprest®

ELISA assay for determination of calprotectin in stools.

For in vitro diagnostic use.

Intended use

Calprest is a quantitative enzyme-linked immunoassay for the detection of calprotectin in faeces. New calprotectin extraction solution is a medium for extraction of calprotectin from faecal samples according to the new and improved faeces sample preparation method [1].

Test Principle

Calprest uses polyclonal antibody against calprotectin in an enzyme linked immuno-sorbent assay system. Calprotectin present in the diluted sample is bound by the antibody adsorbed to the surface of the plastic well. The enzyme conjugated antibody binds to the captured antigen and subsequently the enzyme catalyses the conversion of the substrate to a coloured product. The intensity of the colour is proportional to the amount of conjugate bound, and thus to the amount of captured calprotectin. Concentration of calprotectin in the samples is calculated using the provided standards.

Analytical sensitivity

The analytical sensitivity of "Calprest" is 6.25 ng/ml, which corresponds to 15.6 mg calprotectin/kg faeces at a sample dilution of 1:2500.

Materials provided with the kit

(quantity sufficient for 96 tests)	
Antibody coated plate	12x8 wells
Microtiter plate cover	1 piece
Enzyme conjugate antibody (IgG, rabbit)	1x15 ml
Substrate	1x15 ml
Washing solution (20x)	1x50 ml
Diluent solution(10x)	1x20 ml
Extraction solution(2,5x)	2x50 ml
Standards	6x1 ml
Control 1 (high)	1x 1.0 ml
Control 2 (low)	1x 1.0 ml

Composition of supplied reagents/materials

Antibody coated plate

12x8 wells coated with antibody against calprotectin. Plastic sealed bag containing a desiccant.

Enzyme conjugate antibody (IgG)

1 vials containing 15 ml alkaline phosphatase-labelled anti-human calprotectin IgG antibodies (rabbit) in a buffer solution with the addition of a purple dye. Ready-to-use.

Substrate

1 vial containing 15 ml substrate reagent with Sodium-azide as a preservative. Ready to use.

Washing solution (20x)

1 vial containing 50 ml concentrated washing solution. To be diluted with distilled water.

Dilution solution(10x)

1 vial containing 20 ml diluent solution (10x) with the addition of a blue dye, to be diluted with distilled water.

Extraction solution(2,5x)

2 vials containing 50 ml (2,5x) to be diluted with distilled water.

Standards

6 vials containing 1 ml Calprotectin at six known concentrations (6.25, 12.5, 25, 50, 100 and 200 ng/ml) with the addition of a red dye. The value of each standard is printed on the vial label. Ready to use.

Control 1

1 vial containing 1,0ml with the addition of a red dye. Ready to use. Do not dilute. The range of values is printed on the vial label.

Control 2

1 vial containing 1,0ml with the addition of a red dye. Ready to use. Do not dilute. The range of values is printed on the vial label.

Materials required but not provided

Faeces sample collection

1. Sample collection tube
2. Transport container

Faeces preparation

1. Disposable, breakable sterile inoculation loops
2. Disposable polystyrene screw cap tubes, 14 ml
3. Eppendorf tubes (1 ÷ 1.5 ml)
4. Sensitive digital scale (40 ÷ 150 mg)
5. Vortex mixer
6. Shaker
7. Microcentrifuge (10000g)
8. Freezer (-20°C)

Equipment for ELISA measurements

1. Multi-channel pipette, 50-200 ml
2. ELISA plate washer
3. ELISA plate reader (filter 405 nm)
4. Distilled water
5. Stop solution (NaOH 1M)

Performance criteria

1. Intra-assay

8 reference stool samples were tested 4 times in the same run. The results are reported below:

Samples	1	2	3	4	5	6	7	8
Mean (DO)	0,273	0,509	0,352	0,555	0,201	0,170	0,151	0,744
S.D.	0,011	0,021	0,012	0,015	0,006	0,006	0,004	0,028
CV%	4,0	4,0	3,5	2,7	3,1	3,7	2,7	3,8

2. Inter-assay variation

5 reference stool samples were tested 3 times in 4 different runs. The results are reported below:

Sample	1	2	3	4	5
Mean (ng/ml)	203,00	25,00	38,60	148,70	89,10
S.D.	1,79	1,79	1,44	3,85	3,03
CV%	0,88	7,14	3,74	2,59	3,40

3. Inter-lotto variation

5 reference stool samples were tested in duplicate using 3 different batches. The results are reported below:

Sample	1	2	3	4	5
Mean (ng/ml)	14,13	22,50	12,38	23,61	15,10
S.D.	0,58	1,75	1,05	1,55	0,96
CV%	4,11	7,79	8,47	6,56	6,35

4. Diagnostic sensitivity

Stool from 36 clinically diagnosed IBD pediatric patients were tested. The found value was **95%**.

5. Diagnostic specificity

Stool from 14 subjects with gastrointestinal symptoms not related to inflammatory bowel disease were tested. The found value was **93%**.

Precautions and warnings

1. For in vitro use only.
2. Reagents, samples and microtiter strips should be allowed to reach room temperature (20-25°C) before starting the test.
3. Materials of human origin used in this kit have been tested and confirmed negative for HBsAg and anti-HIV I and II and anti-HCV antibodies. However, they should be treated as a potential biohazard, and shall be handled and disposed of according to local laboratory legislation.
4. **Warning:** do not interchange components from the different kit batches. Satisfactory performance of the test is guaranteed only when components from the same batch of Calprest are used.
5. The substrate reagent contains sodium azide as preservative at concentrations lower than 0,1% (w/w). Sodium azide is highly toxic to aquatic organisms. Avoid to empty into drains.
6. Substrate reagent should be pale yellow. The substrate is light sensitive. Store in the dark and shake before use.
7. Unused microtiter strips should be re-sealed air tight in the foil bag with the enclosed drying pad and stored at 2-8°C.
8. Insufficient washing of the ELISA plate can lead to erroneous values of Calprotectin due to incomplete removal of reagents. Routine maintenance of aspiration/wash system is strongly recommended.
9. All reagents, except to the substrate and the concentrated washing solution, contain Proclin 300 as a preservative agent below the allowed limits.

Reagent preparation

Extraction buffer

Dilute concentrated Extraction solution by adding 1 part (50 ml) of it to 1,5 parts (75 ml) of freshly distilled water to obtain 125 ml working solution. Mix well.

Diluent buffer

Dilute concentrated Dilution solution by adding 1 part (20 ml) of it to 9 parts (180 ml) of distilled water to obtain 200 ml working solution. Mix thoroughly.

Washing solution

Prepare the washing solution by diluting the content of the whole vial (50 ml) with distilled water to a final volume of 1000 ml.

Storage

All reagents must be stored at 2-8°C and used before expiration date printed on the label.

Reagent stability (unsealed reagents)

Reagent	Storage conditions	Storage time
Conjugate	2-8°C	1 month
Substrate	2-8°C	3 months
Standards	2-8°C	1 mese
Controls	2-8°C	1 mese

Stability of working solutions

Reagent	Storage conditions	Storage time
Washing solution	20-25°C	30 days
Extraction solution	2-8°C	3 months
Dilution liquid	2-8°C	30 days

1. The expiry date is printed on all component labels.
2. Avoid exposure to high temperature, direct sunlight or extreme humidity.
3. Unused microtiter strips resealed air tight in the plastic bag with the drying pad inside can be stored for one month.

Stability during transportation

An accelerate stability test proved that all reagents are stable after 96 hours at 37°C.

Faeces sample collection and preparation

1. Collect approx. 1 ÷ 5 g stool and place it in a suitable container.
2. Stool sample is sent or delivered to the laboratory within 4 days.
3. Temperatures during transport should never exceed 30°C.
4. Store faeces samples at -20°C.
5. Thaw frozen stool samples at room temperature.
6. Weigh (tare) the empty screw cap tube together with the inoculation loop.
7. Take out approx. 100 mg (between 40 - 120 mg) faeces by means of the inoculation loop, and place into a screw-cap-tube.

8. Weigh tube and loop with faeces and calculate net faeces weight (between 40 - 120 mg).
9. Break off the loop handle, leaving the loop with faeces and a 4-6 cm handle inside the screw cap tube.
10. Add pre-diluted extraction buffer (weight/volume ratio 1:50), e.g. 100 mg faeces + 4.9 ml diluted extraction solution. (See table in the Technical advice and customer service). Close the tube.
11. Shake/mix vigorously for 30 seconds by means of a vortex.
12. Homogenise 25 ± 5 minutes on a shaker or roller. The loop inside the tube will act as an agitator.
13. Transfer the homogenate (1 ml) to an Eppendorf tube and centrifuge for 20 min. at 10,000g at RT using a bench-top centrifuge (e.g., Heraeus Biofuge 13).
14. Transfer 0.5 ml of the clear extract supernatant to new Eppendorf tubes.
15. The extracts may be tested immediately. Store frozen extracts (max. 3 months at -20°C) for later measurement.

ELISA Procedure

The procedure described hereafter can be used for testing plasma samples. See "Technical Advice" for further details.

1. Ensure that all reagents reach room temperature (20-25°C).
2. Thaw frozen sample at room temperature.
3. Dilute samples 1:50 (20 µl sample + 980 µl dilution buffer). For further dilution of high concentration samples, dilute the sample to a final 1: 250 dilution (e.g. 200 µl of the 1:50 dilution + 800 µl dilution buffer).
4. Suggested plate layout in duplicates is shown below. Fit the strip holder with the required number of micro ELISA strips. Use uncoated strips to complete the strip holder if the washer requires a full plate. Blank, standards and controls must be included in each run.

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5. Add 100 µl of Dilution buffer to wells A1-B1 (blank).
6. Add 100 µl of each standard in duplicate wells (C1-D1, E1-F1, G1-H1, A2-B2, C2-D2, E2-F2).
7. Add 100 µl of each control in duplicate wells (G2-H2, A3-B3).
8. Mix diluted sample well before application to the plate and add 100 µl of each sample in duplicate wells (C3-D3, E3-F3...).
9. Cover plate with plate cover and incubate the plate at room temperature for 45 ± 5 min.

10. At the end of the incubation time, wash the plate by adding 0,3 ml of diluted washing solution to each well. Remove as much liquid as possible. Repeat this step two more times up to a total of 3 washing steps. After the final aspiration, invert plate and tap gently on absorbent tissue to ensure complete removal of washing solution.
11. Add 100 µl conjugate to each well.
12. Cover plate with plate cover and incubate the plate at room temperature for 45 ± 5 min.
13. Repeat washing step as above (see 10).
14. Add 100 µl substrate solution to each well. Note: a multi-channel pipette is recommended in order to avoid variation in substrate development time. Avoid also the formation of bubbles during substrate pipetting.
15. Incubate the plate at room temperature for approx. 30 minutes in a dark place or wrap the plate with an aluminium foil.
16. Read the O.D. values by means of an ELISA reader at 405 nm. When standard 6 reaches an OD value between 1.8-2.1, the reaction should be read with an automatic EIA reader or stopped by adding 100 µl NaOH 1M stop solution. Plates stopped with NaOH 1M may be stored at 4°C for 24 hours.

Assay calibration and quality control

1. A new standard curve is used with each run.
2. The plate background should be < 0.20 OD.
3. Control 1 and 2 are to be included in each run.

Calculation of test results

Calculate the mean optical densities (OD) of all duplicates. Subtract the mean Blank OD from all values. Plot the standard curve with actual calprotectin concentration of standards as ng/ml and corresponding mean OD values on an xy system. The readings of the Samples from the Standard curve is corrected for the dilution and converted to mg/kg by multiplying 2.5 (e.g. a reading of 50 ng/ml becomes 125 mg/kg). If samples are further diluted this must be compensated for during calculations. Concentrations can also be determined by use of a computer linked to the Elisa reader.

Reference values

Clinical studies^(2-5, 17) gave the following values:

Normal healthy reference, median	25 mg/kg
Colorectal cancers, median	350 mg/kg
Inflammatory bowel disease (CD&UC), median	1722 mg/kg

Interpretation of results

Samples giving values above 50 mg/kg are regarded as having a positive Calprest test. In healthy adults the median value is about 25 mg/kg. The median value in patients with colorectal cancers is about 350 mg/kg. Patients with active, symptomatic inflammatory bowel disease have levels between 200 and 20 000 mg/kg.

References

1. Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
2. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992;27:793-798.
3. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
4. Tibble J. et al.: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000;47:506-513.
5. Bunn S.K. et al.: Fecal calprotectin: Validation as a noninvasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001;33:14-22.
6. Dale I. et al.: Purification and partial characterization of a highly immunogenic human leukocyte protein, the L1 antigen. Eur J Biochem 1983;134:1-6.
7. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.
8. Brandtzaeg P. et al.: Distribution of a formalin resistant myelomonocytic antigen (L1) in human tissues. II. Normal and aberrant occurrence in various epithelia. American J of Clin Pathology 1987;87:700-707.
9. Fagerhol M.K.: Nomenclature for proteins: is calprotectin a proper name for the elusive myelomonocytic protein? J Clin Pathol: Mol Pathol 1996;49:M74-M79.
10. Isaksen B. and Fagerhol M.K.: Calprotectin inhibits matrix metalloproteinases by sequestration of zinc. J Clin Pathol: Mol Pathol 2001;54:289-292.

11. Steinbakk M. Et al.: Antimicrobial actions of calcium binding leukocyte L1 protein, calprotectin. Lancet 1990;336:763-765.
12. Yui S. et al.: Induction of apoptotic cell death in mouse lymphoma and human leukemia cell lines by a calcium-binding protein complex, calprotectin, derived from inflammatory peritoneal exudate cells. Journal of Leukocyte Biology; 58:650-658
13. Johne B. et al.: Functional and clinical aspects of the myelomonocyte protein calprotectin. J Clin Pathol: Molec Pathol 1997;50:113-123.
14. Ton H. et al.: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000;292:41-54.
15. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000;95:2831-2837
16. Røseth A.G. et al.: Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.
17. Johne B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

Calprest

96 tests, code 9031

Date of preparation:

2012.10.25

Eurospital

Eurospital SpA

34147 Trieste, via Flavia 122

Tel. +39 040 8997.1 Fax +39 040 280944

www.eurospital.com

Sample preparation

a) Stools

Weighing and sample dilution

For a rapid sample dilution regime the table below gives the volume of extraction solution added to a given amount of stools (weight + volume = 1+49)-

Stools (mg)	Extraction solution (ml)
120	5.9
115	5.6
110	5.4
105	5.2
100	4.9
95	4.7
90	4.4
85	4.2
80	3.9
75	3.7
70	3.4
65	3.2
60	2.9
55	2.7
50	2.5
45	2.2
40	2.0

Extraction and homogenisation

1. Vortex mixing time should be as accurate as possible (30 seconds).
2. Shaking time not shorter than 20 minutes. Reduce delay before pipetting.
3. Recommended batch sizes: 10, 20, 30 or 40 samples
4. Centrifugation: for ease of harvesting and handling, make sure the pellet is compact at the bottom of the centrifuge tube.
5. Do not change centrifugation time and g force. Check g force and rpm.

Harvesting of faeces extract supernatant

1. Harvest the top half of the clear supernatant, approx. 500 µl.
2. Avoid contact with the pellet as aggregates or particles can cause erroneous calprotectin values.

b) Plasma

The ELISA procedure is the same for both faeces and plasma specimens.

Materials required but not provided

1. EDTA Vacutainer® Tubes for plasma collection
2. Centrifuge
3. Eppendorf tubes

Sample collection

1. Draw blood into EDTA Vacutainer® tube
2. Centrifuge 2000g for 15 minutes
3. Plasma is pipetted off and frozen at -20°C in aliquots larger than 150 ml

Calculation of test results

The calprotectin concentration in the sample is determined by multiplying the actual calprotectin concentration read from the standard curve by 0,05. If samples are further diluted, this factor must be included in the calculation. The concentration is expressed as µg/ml plasma.

Reference values

A preliminary study gave the following values:

Normal reference range: 0,3-1,6 µg/ml

Note on ELISA assay procedure

Washing

1. Avoid blocking of aspiration probes

Pipetting

1. Use disposable tips
2. Avoid contamination of the pipette
3. Avoid air bubbles in the tips

Substrate

1. Store Substrate container in the dark.
2. Colour of substrate should be colourless or pale, not yellow.
3. Use separate reservoir for conjugate and substrate if different than original vials.
4. Change vessel between each plate. Use multichannel pipettes to eliminate variations in substrate development time. Carry-over from conjugate to substrate is avoided if only one pipette is used.

Conjugate

1. Mix content of the vial gently prior to use (do not shake).

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Références

- Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
- Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992;27:793-798.
- Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
- Tibble J. et al.: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000;47:506-513.
- Bunn S.K. et al.: Fecal calprotectin: Validation as a noninvasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001;33:14-22.
- Dale I. et al.: Purification and partial characterization of a highly immunogenic human leukocyte protein, the L1 antigen. Eur J Biochem 1983;134:1-6.
- Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.
- Brandtzaeg P. et al.: Distribution of a formalin resistant myelomonocytic antigen (L1) in human tissues. II. Normal and aberrant occurrence in various epithelia. American J of Clin Pathology 1987;87:700-707.

9. Fagerhol M.K.: Nomenclature for proteins: is calprotectin a proper name for the elusive myelomonocytic protein? J Clin Pathol: Mol Pathos 1996;49:M74-M79.

10. Isaksen B. and Fagerhol M.K.: Calprotectin inhibits matrix metalloproteinases by sequestration of zinc. J Clin Pathol: Mol Pathol 2001;54:289-292.

11. Steinbakk M. Et al.: Antimicrobial actions of calcium binding leukocyte L1 protein, calprotectin. Lancet 1990;336:763-765.

12. Yui S. et al.: Induction of apoptotic cell death in mouse lymphoma and human leukemia cell lines by a calcium-binding protein complex, calprotectin, derived from inflammatory peritoneal exudate cells. Journal of Leukocyte Biology; 58:650-658

13. Johnne B. et al.: Functional and clinical aspects of the myelomonocyte protein calprotectin. J Clin Pathol: Molec Pathol 1997;50:113-123.

14. Ton H. et al.: B. Improved assay for fecal calprotectin. Clinica Chimica Acta 2000;292:41-54.

15. Limburg P.J. et al.: Fecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000;95:2831-2837

16. Røseth A.G. et al.: Assessment of disease activity in ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.

17. Johnne B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

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8. Brandtzaeg PÅG. y otros: Distribution of a formalinresistant myelomonocytic antigen (L1) in human tissues. II. Normal and aberrant occurrence in various epithelia. American J of Clin Pathology 1987;87:700-707.

9. Fagerhol M.K.: Nomenclature for proteins: is calprotectin a proper name for the elusive myelomonocytic protein? J Clin Pathol: Mol Pathos 1996;49:M74-M79.

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12. Yui S. y otros: Induction of apoptotic cell death en mouse lymphoma and human leukemia cell lines by a calcium-binding protein complex, calprotectin, derived from inflammatory peritoneal exudate cells. Journal of Leukocyte Biology; 58:650-658

13. Johnne B. y otros: Functional and clinical aspects of the myelomonocytic protein calprotectin. J Clin Pathol: Molec Pathol 1997;50:113-123.

14. Ton H. y otros: B. Improved assay for fecal calprotectin. Clinica Chimica ACTA 2000;292:41-54.

15. Limburg PÅG.J. y otros: Faecal calprotectin levels predict colorectal inflammation among patients with chronic diarrhea referred for colonoscopy. Am J Gastroenterol 2000;95:2831-2837

16. Røseth A.G. y otros: Assessment of disease activity en ulcerative colitis by faecal calprotectin, a novel granulocyte marker protein. Digestion 1997;58:176-180.

17. Johnne B. y otros: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

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Literatur

1. Fagerhol M.K. et al.: Calprotectin (The L1 leukocyte protein) in: Smith V.L. and Dedman J.R. (eds.): Stimulus response coupling: The role of intracellular calcium-binding proteins. CRC Press, Boca Raton 1990, p.187-210
2. Røseth A.G. et al.: Assessment of the neutrophil dominating protein calprotectin in feces. Scand J Gastroenterol 1992;27:793-798.
3. Røseth A.G. et al.: Correlation between faecal excretion of Indium-111-labelled granulocytes and calprotectin, a granulocyte marker protein, in patients with inflammatory bowel disease. Scand J Gastroenterol 1999; 34:50-54.
4. Tibble J. et al.: A simple method for assessing intestinal inflammation in Crohn's disease. Gut 2000;47:506-513.
5. Bunn S.K. et al.: Fecal calprotectin: Validation as a noninvasive measure of bowel inflammation in childhood inflammatory bowel disease. J Pediatr Gastroenterol Nutr 2001;33:14-22.
6. Dale I. et al.: Purification and partial characterization of a highly immunogenic human leukocyte protein, the L1 antigen. Eur J Biochem 1983;134:1-6.

7. Sohnle P.G. et al.: The zinc-reversible antimicrobial activity of neutrophil lysates and abscess fluid supernatants. Journal of Infectious Diseases 1991;164:137-142.

8. Brandtzaeg P. et al.: Distribution of a formalin resistant myelomonocytic antigen (L1) in human tissues. II. Normal and aberrant occurrence in various epithelia. American J of Clin Pathology 1987;87:700-707.

9. Fagerhol M.K.: Nomenclature for proteins: is calprotectin a proper name for the elusive myelomonocytic protein? J Clin Pathol: Mol Pathos 1996;49:M74-M79.

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17. John B. et al.: A new fecal calprotectin test for colorectal neoplasia. Clinical results and comparison with previous method. Scand J Gastroenterol 2001;36:291-296.

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3. Reaktionsgefäße, z.B. der Firma Eppendorf AG